

## TROUBLES AND TROUBLESHOOTING

<b>Troubles</b>	<b>Checklists</b>	<b>Countermeasures</b>
The measurement is not accurate	Please check if the infrared sensor is clean or not	Clean the infrared sensor with tissue or dry soft cloth gently
	Please check if the way how the measurement is done is correct or not	Ensure you have read the manual and know how to use the product properly.
	Please check if you have let the thermometer gets stabilized in the room for 30 minutes at least	Please keep the thermometer in the room for at least 30 minutes before usage (precision : $\pm 1^{\circ}\text{C}$ )
	Are you using the thermometer indoor ?	Please take the measurement indoor
	Please check if you have hold the thermometer in your hand too long and affect the accuracy	Put the thermometer on the table in the room where the measurement is taking place and let it cool down first (precision : $\pm 1^{\circ}\text{C}$ )
	Please check if there is low battery icon on LCD	Replace with new batteries